

September 5, 2017

SHELF LIFE FOR ECG ELECTRODES IN OPEN POUCH

Dear Valued Customer:

ECG-electrodes manufactured by TYROLMED GmbH feature either a wet gel which is designed to stay fresh for up to 14 days or a solid gel which is designed to stay fresh for up to 30 days when stored in an opened bag or tray, provided storage conditions given on the package are followed.

For additional technical questions, please contact our professional customer service or sales representative.

Sincerely,



Ing. Mag. Lukas HAFFNER
Quality & Regulatory Affairs Manager
Tyrolmed GmbH, Austria