




FEATURES

- 5-day wear time means fewer disruptions ensuring the body is focused on essential growth and repair
- Lightweight leadwires for use on even the tiniest of patients
- Pre-wired electrodes helps to eliminate motion artifact

Lightweight, disposable, pre-wired electrodes for the NICU | PICU



PACKAGING

Electrodes are packaged in a single-use pouch with an easy-to-open tear

Single-Use Pouch




FRONT



BACK

Box





MONITORING

Product #	Ref #	Description	Gel Type	Duration of Use
4009832C	50401 - N3	kit of 3 electrodes for neonates		5 Days

COMPATIBLE MONITOR CONNECTIONS

Compatible Cables include adaptors for Bard, Biosense, Mindray, Datascope, GE Marquette, Philips, Datex, Mennen, Mortara, Nihon Kohden, GE, Getemed, Spacelabs, Siemens, Drager, Schiller, Zoll and more.



Mindray



Philips



Spacelabs
TruLink



GE



Physio



Nihon
Kohden

3-13-01-049 (06/2017)

CABLE CONNECTIONS

Electrode → Patented Electrode Cable Connector → Cable → Monitor

