



FEATURES

Our high quality ultrasound and ECG gels are specially formulated to ensure excellent signal transmission.

- Water soluble
- Odourless
- Non-irritant
- Oxidants free
- Non-greasy
- Salt free
- Alcohol free
- Abrasive free



PACKAGING



250 ml Bottle
12 / Case

5 litre Bottle
4 / Case
with refillable
dispenser and cap

+44 1803 860100 | Vermed.co.uk

Patient Monitoring & Medical Surgical Consumables

Vermed Ultrasound Gels

Standardization Made Simple™

VERMED ULTRASOUND GEL

Product #	Description	# / Box
4009235C	250 ml Bottle - Blue Ultrasound Gel	12
4009442C	5 litre Bottle - Blue Ultrasound Gel	1
4009507C	250 ml Bottle - Clear Ultrasound Gel	12
4009508C	5 litre Bottle - Clear Ultrasound Gel	1



WE ALSO CARRY

Thermasonic Gel Warmers

Single bottle and three bottle warmer heats ultrasound gel to body temperature. Durable wall mountable construction, with warranty.

- Powerful heating element uniformly heats bottles of gel
- Microprocessor control maintains constant temperature
- "Thermal cutoff" provides protection from overheating
- IEC320/C14 receptacle
- Hospital-grade plug and power cord
- High-density, impact-resistant composite
- Adjustable temperature settings for multi-bottle units
- Illuminated power ON/OFF switch on multi-bottle units
- Heat indicator lamp on single bottle unit
- Low electrical demand

Product #	Cross Ref.	Description
4007906C	82-01-20 CE	Single Bottle Warmer - 230V (Euro Plug)
4007908C	82-03-20 CE	Multi Bottle Warmer - 230V (Euro Plug)
4007909C (Pictured Right)	83-03-20 CE	Multi Bottle Warmer LCD - 230V (Euro Plug)

